

HAIRMAX[®] LaserBand 82

QUICK START GUIDE



Use your HairMax for as little as 90 seconds, 3 times a week, ideally every other day (for example, Monday, Wednesday, Friday).

FOR MAXIMUM BENEFIT:

- Scalp should be freshly washed (hair can be wet, dry or damp)
- Use your HairMax prior to applying styling products (these can build up and impede the laser light)
- Regularly clean the laser area with a damp cloth. Do not use harsh detergents or cleaning products on any part of your HairMax

CHARGING

Charge your HairMax for 4-5 hours before initial use.

Prior to charging your HairMax, make sure your HairMax is turned off.




Charge your HairMax by plugging the power cord directly into the charging socket on the bottom right side of your HairMax. See the User Manual for the correct orientation of the power cord (Figure A). Plug the other end into an electrical socket.

The indicator light above the power button will flash during the charging process. A steady light indicates your HairMax is fully charged (Figure B).

Do not use your HairMax while it is plugged into the power cord.

Prior to treatment, test your HairMax for sensitivity. For more information, see the user manual.



Next, press the power button once to activate the 30-second treatment method (recommended) OR twice for the 4-second treatment method (Figure C).



TREATMENT METHOD

30-Second Treatment Method (Recommended)

Press the power button once to select the 30-second treatment method. Your HairMax will vibrate every 30 seconds. Place your HairMax at your hairline and leave it in place until it vibrates (Figure D). Next, slide your HairMax to the middle portion of your scalp using light pressure so the teeth can separate your hair (Figure E). When your HairMax vibrates again, slide it to the last spot on your scalp (Figure F). You have completed today's treatment.




4-Second Treatment Method

Press the power button twice to select the 4-second treatment method. Your HairMax will vibrate every 4 seconds. Place your HairMax at your hairline and leave it in place until it vibrates (Figure D). Next, slide your HairMax to the middle portion of your scalp using light pressure so the teeth can separate your hair (Figure E). When your HairMax vibrates again, slide it to the last spot on your scalp (Figure F). This completes the first pass of your treatment. Repeat this process 6 times to complete today's treatment.

After treatment, turn off your HairMax by pressing and holding the power button until the indicator light turns off (Figure G). Remove your HairMax from your head and plug it into the charger for future use.

NOTE: The LaserBand 82 turns off automatically after 10 minutes.



The HairMax® LaserBand is a highly effective device when used correctly. It is manufactured using precision components and should last many years of normal use. If you experience difficulties with your device, please refer to the troubleshooting section of your user manual.

For detailed information about proper use, warnings and cautions, please take time to read the entire user manual before operating your HairMax.

HairMax® offers alternate teeth sizes (shorter or longer) to maximize the fit and comfort of your device. See the User Manual for more information.

Laser Light - Avoid Direct Eye Exposure

HairMax® laser devices are indicated to treat Androgenetic Alopecia, and promote hair growth in males who have Norwood-Hamilton Classifications of IIa to V and in females who have Ludwig (Savin) I-4, II-1, II-2, or frontal patterns of hair loss and who both have Fitzpatrick Skin Types I to IV.